



# Massager acupuncture mat

9 Vibrating Motors Full Body Far Infrared Heat  
heathy Massage Mat cushion seat mattress

---

## INSTRUCTIONS GUIDE

---

**MODEL BYM620C**

PROVIDED BY WAL-MART CANADA, MISSISSAUGA, ONTARIO L5N 1P9  
PRINTED IN CHINA  
PATENT PENDING

Two sides usable with different modes:



## HOW TO USE THE REMOTE CONTROL?

The massage mat is equipped with an easy-to-use remote control, to use the 9 massagers in 4 body zones plus soothing warmth.

1. Each massage zone is individually powered. To activate a specific massage zone, slide the button of the desired zone to the HI (High) or LO (Low) position.
2. You can activate a massage zone or a combination of zones simultaneously. The LED lights indicate the selected zones.



3. Body massage areas:  
 Zone 1: neck & shoulders  
 Zone 2: Upper & Lower Back  
 Zone 3: Legs  
 Zone 4: Calf
4. Slide all massage buttons to the OFF position before removing the AC plug from the wall outlet.
5. To turn on the heat, press the red heat button. To deactivate, press the heat button again. The LED indicator light will illuminate when the heat mode is ON.
6. To disconnect, set the control buttons to OFF, then remove the AC adapter from the wall outlet.

## AC WALL ADAPTER

1. To turn on the massage mat, plug the power cord into the lower part of the controller.
2. Make sure that all control buttons are in the OFF position before plugging the AC adapter into the wall outlet.
3. To activate, slide the control keys to the HI (Up) to LO (Down) position.
4. To disconnect, slide the control buttons to OFF to remove the AC adapter from the wall outlet.



## ELASTIC STRAPS

1. Place the massage mat on your favorite chair and secure it securely with the handy elastic straps. You can also attach a cushion with the elastic straps to raise your neck and shoulders for a relaxing massage.
2. The reversible faux-sheepskin mat can be rolled up using the convenient elastic straps as shown. This roll-up design is easy to wear and store.

## IMPORTANT SAFETY INSTRUCTIONS

(PLEASE KEEP IT A SAFE PLACE)

### DANGER TO REDUCE THE RISK OF ELECTRIC SHOCK:

1. Always unplug this appliance from the wall outlet immediately after use or before cleaning.
2. DO NOT bend over to pick up an appliance that has fallen into water. Unplug immediately.
3. DO NOT use while taking a bath or shower.
4. DO NOT place or store the appliance where it could fall or end up in a bathtub or sink.
5. DO NOT place or put in water or other liquids.

### WARNING

TO REDUCE THE RISK OF BURNS, FIRE, ELECTRIC SHOCK OR PERSONAL INJURY:

1. This appliance should not be left unattended when plugged in. Unplug it from the wall outlet.
2. DO NOT use under a blanket or cushion. It may overheat and cause fire, electric shock, or injury to persons.
3. Special supervision is required when the appliance is used by children, handicapped or disabled persons. THIS APPLIANCE IS NOT A TOY.
4. Use only this unit for its intended use as described in this manual.
5. DO NOT use this appliance if it has a damaged power cord or plug, if it is not working properly, if it has fallen to the ground and sustained damage, or if it has been dropped into water.

6. DO NOT carry this appliance by the power cord or use the cord as a handle.
7. Place the appliance and cord away from heated surfaces.
8. DO NOT USE outdoors .
9. Never place or insert any object into any opening of the appliance.
10. DO NOT use in areas where aerosol products (sprays) are used.
11. To disconnect, turn all buttons to the OFF position, then unplug the unit.
12. Never use while sleeping and avoid falling asleep while the product is running.
13. Use only the original adapter that was sold with the unit.

## **ATTENTION**

THE DEVICE SHOULD NOT BE USED IN THE FOLLOWING CASES:

1. If you have any health concerns, please consult your physician before using this product.
2. Persons with heart pacemakers should consult a physician before using this product.
3. Never leave the massage mat unattended, especially in the presence of children.
4. DO NOT fall asleep while using this massager.
5. The massage should be pleasant and enjoyable. If you experience pain or discomfort after use, stop using it and consult your physician.
6. DO NOT place the remote control under the seat.
7. Avoid using an extension cord, which could overheat and cause a fire hazard.
8. Do not jump or fall suddenly on this product.
9. Do not use this product before lying down. Massage could create a stimulating effect and delay sleep.

## **CARE AND MAINTENANCE:**

1. No maintenance of this product is required. There are no spare parts available for the user.
2. Store the seat set in a cool, dry place, away from moisture when not in use.
3. DO NOT wrap the power cord around the unit.
4. DO NOT hang the unit by its power cord.
5. Keep any strong solvents or detergents away.
6. Unplug the unit and allow to cool before cleaning. Wipe with a soft, slightly damp cloth. Never immerse the unit in water or other liquids.
7. Do not use abrasive cleaners, brushes, gasoline, etc. to clean the massage mat.
8. Regularly brush the faux sheepskin with a soft brush and vacuum occasionally.
9. If the carpet is wet, simply air dry.
10. DO NOT wash in the washing machine.