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VIBRO PLATE INSTRUCTION

TROUBLE SHOOTING	-----	2
IMPORTANT SAFETY INFORMATION	-----	2-3
PRODUCT DISCRIPTION	-----	4
DISPLAY INSTRUCTION	-----	5
EXERCISE SESSION	-----	6

TROUBLE SHOOTING & IMPORTANT SAFETY INFORMATION

Please keep this manual in a safe place for references.

WARNING: Please read carefully the following important precautions and information before using the unit for the first time to avoid the risk of burns, electric shock, fire or injury.

Trouble shooting

1. Emergency stop: LED on Console shows “ER-3”, LED flashes three times on PCB.
2. Connection error: LED on console shows “ER-4”, LED flashes four times on PCB.
3. Electron Power tube damage: LED on console shows “ER-5”, LED flashes five times on PCB.
4. Current Overload: LED on console shows “ER-6”. LED flashes six times on PCB.
5. Engineer abnormality: LED on Console shows “ER-7”, LED flashes seven times on PCB.
6. Blocked Engineer: LED on Console shows “ER-8”, LED flashes eight times on PCB.
7. Low voltage: LED on Console shows “ER-9”, LED flashes nine times on PCB

Import safety information

1. It is the responsibility of the owner to ensure that all users of this device are adequately informed about the operation and their precautions.
2. Use this device only as described in the manual
3. Always keep the device dry and dust-proof.
4. Do not operate the device while using aerosol or administering oxygen
5. Place the device on a sufficiently level surface (at least 8 feet). It is recommended to put a mat under it.
6. The device is only suitable for one person at the same time.
7. Keep children under 12 years and pets away from the unit.
8. Only suitable for persons up to 120 KG
9. Keep the power cord and surge protector away from hot surfaces.
10. Check all parts of the device at regular intervals.

11. Never leave the device unattended during startup and remove all cables when not in use.
12. Do not try to adjust or move the unit unless it is mounted properly.
13. Never insert or drop an object into an opening of the unit.
14. **DANGER:** Disconnect the power cord after use before performing any of the maintenance and adjustment procedures described in the manual. Only open the bonnet when instructed to do so by an authorized service representative.
15. Maintenance that is not in the manual is only carried out by authorized customer service representatives.
16. This device is intended for home use only. Use in commercial, leased or institutional facilities is prohibited.

WARNING: It is strongly recommended that you consult your doctor before first use. This applies in particular to persons with health problems as well as to persons over 35 years of age. Read all instructions before first use

DISPLAY DESCRIPTION

PRODUCT DATA

Size	67.0 x 38.0 x1 2.5 CM
Weight	12.5 KGS
Power	200W
Input	220V 50Hz
Frequency	5 to 15Hz
Control mode	PWM CONTROLLER [DC MOTOR DRIVER]

COMPONENTS

PANEL DISPLAY INSTRUCTION

Power: Power on or off

Prg Mode: Select the Program model

H.R: Manual Model

You can adjust the time and speed separately before or after machine works

P0-P9: Automatic Model


- 1.You can adjust the time before or after working.But Speed can not do this.
- 2.The speed will automatically vary as a certain rhythm.
- 3.Different automatic model,the rhythm is different.

Time: This can increase or reduce time(0-15Minutes)

Speed: This can increase or reduce speed(0-99levels).

Start/Stop: Press this, the machine begin/stop to work.

REMOTE CONTROLLER INSTRUCTION

EXERCISE SESSION

The fitness vibration trainer uses the body's own reflexes to exercise your muscles in a convenient manner that does not require that you get out of breath. The fitness vibration trainer accomplishes this best when you are in a stress position, such as a squat posture, and your muscles are working to keep you in that position. When the oscillating plate moves quickly, your muscles automatically move a small amount to keep you in the posture, hence exercising your muscles really efficiently. You will find the fitness vibration trainer extremely effective when using it in one of the positions shown below. You will get the best results if you bend your knees or arms (depending which is pressing on the oscillating plate) so your muscles adopt the automatic motion described.

You can also use the fitness vibration trainer as a general massager to ease muscle pain and improve circulation. Here, you should not press so hard on the plate, rather let the plate do the work for you. Use several different postures for a few minutes each for a complete workout.

			
<p>Whole-Body-Exercise <i>Stand on the plate and stretch your feet as wide as your shoulders. This will help improve overall circulation and muscle tone.</i></p>	<p>Squat Posture <i>Stretch your legs as wide as your shoulders when squatting and bend your knees down to 90 degrees, if you can. Hold the position. This is great for a thigh workout.</i></p>	<p>Press-Up Posture <i>Stretch your arms as wide as your shoulders. If you wish to exercise the muscles more, bend your elbows to 90 degrees and hold this position.</i></p>	<p>Lower Legs Posture <i>Put your lower legs as shown on the plate with your hand supporting behind you on the ground. For an even more effective exercise, raise your seat off the floor and press down on your calves.</i></p>
			
<p>Waist Bend Posture <i>Stretch your legs, bend at the waist and put your hands on the plate. Bend your elbows for a more effective exercise.</i></p>	<p>Single-Foot Posture <i>Put one foot on the pedal and relax yourself. For a more effective exercise, lunge forward, placing more body weight on the plate.</i></p>	<p>Sitting Posture 1 <i>Sit with your buttocks on the oscillating plate. Sit comfortably, keeping your back straight.</i></p>	<p>Sitting Posture 2 <i>Put your legs on the plate when sitting on a chair. This is a great way to ease tension and can help to improve circulation in the legs.</i></p>

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